


## US Endovascular .014/.018 Nitinol Interventional Guide Wire

USE GS is an innovative .014/.018 nitinol interventional guide gold wire for successful and quick procedures in lower limb ischemia

### Highlights

- Nitinol kink-resistant core with a radiopaque Galt Gold® tip for enhanced visibility.
- Ideal for complex cases, facilitating smooth exchanges and navigating tortuous vessels without losing lubricity during high friction procedures.
- Hydrophobic silicon coating compatible with atherectomy devices.
- 1:1 torque and steerability.
- Engineered with transition-less nitinol for smooth delivery.
- 5cm radiopaque gold tip for enhanced performance and visibility.
- Available in either .014/.018 180cm/300cm lengths.



### FDA 510k Cleared Sizes

Guidewire Length (cm)	All Sizes Available in Straight or Angled Tip			
	.014 Guide Wire Diameter	.018 Guidewire Diameter	Straight Tip	Angled Tip
180	✓	✓	✓	✓
300	✓	✓	✓	✓